
合作创造美好未来  
Collaboration Makes a Better Future

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## Hollow Fiber Cartridges User Manual

Please contact us for more information  
Hangzhou Cobetter Filtration Equipment Co.,Ltd.

Sales Add: 13-18F, Building 1, Zicheng International Innovation Center,  
No.39 Jincheng Road, Xiaoshan District, Hangzhou 311215, China  
Factory Add: Cobetter Park, Heshang New Material Industry Park,  
Xiaoshan District, Hangzhou 311625, China

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使用前请仔细阅读  
Please read the user manual carfully.


## 1. Preparation of Hollow Fiber

For laboratory or pilot scale application, the basic manual control system includes a pump, a feed tank, a permeate collection tank, pressure gauges and valves, and the system is connected with hoses and connectors of appropriate specifications.

The pressure gauges shall be filled with glycerin or mechanically wetted out.

Usually, a hose clamp is installed at one permeate port to control the outflow of the permeate, and the other permeate port is set closed. A pressure gauge may be connected to the permeate end if necessary.

Cobetter's hollow fiber module initially contains glycerin to protect the membrane structure, which needs to be rinsed before use. First, add sufficient deionized(DI) water to the feed tank, close the permeate valve, open the retentate valve, discharge waste via the retentate line, and maintain an inlet pressure of 2-5psi. The volume of flushing liquid should be greater than 20L/m<sup>2</sup>. Then, move the retentate line back to the feed tank, discharge the waste through the permeate end, open the permeate valve, adjust TMP to 2-5psi, and the permeate volume should also be greater than 20L/m<sup>2</sup>.

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

## 2. Sterilization

If necessary, sanitize after the system has been thoroughly cleaned. The recommended disinfectant is 0.1M NaOH solution.


1. Add 10-20L/m<sup>2</sup> disinfectant to the feed tank. Open the retentate valve and the permeate valve, and place the retentate end and the permeate end into the feed tank.
2. Adjust the pump speed and retentate valve to maintain a TMP of 2-5psi and circulate for 30min.
3. Drain the disinfectant in the system after circulation.
4. Add sufficient DI water to the cleaned feed tank, close the permeate valve, open the retentate valve, direct the retentate line to drain, adjust the pump speed to maintain the inlet pressure at 2-5psi, and flush until the outflow liquid is neutral.
5. Move the retentate line into the feed tank, open the permeate valve and direct the permeate line to drain, adjust the pump speed and the retentate valve to maintain TMP at 2-5psi, and process until conductivity or pH of the flushing liquid meets customer's requirements.

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
### 3. Measuring Water Flux

1. Add enough DI water to the feed tank, open the permeate valve and retentate valve, place both the permeate end and retentate end into the feed tank, and drain the gas in the system.
2. Adjust the pump speed and retentate valve to maintain the TMP at 3 psi, then measure and record the permeate flow rate and water temperature in the inlet tank.
3. Adjust the pump speed and retentate valve to maintain the TMP at 7 psi, again measure and record the permeate flow rate and water temperature in the inlet tank.

4. Use the following equation to calculate water flux.

$$\text{Flux in LMH} = \frac{\text{Permeate flow (mL/min)} \cdot 0.06}{\text{Cassette surface (m}^2\text{)}} \quad (\text{l/m}^2/\text{hr})$$

5. Multiply the calculated flux by temperature correction factor for corresponding temperature from Table 1 to obtain standard water flux, and then divide by corresponding pressure value to obtain the unit pressure water flux (per psi).
6. If there is no large deviation between the two values, the average of the two may be taken as the standard water flux of the module.
7. If there is a large deviation between the two values, repeat steps 2-6.

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
Table 1. Viscosity correction factor

Temperature in °C when permeate flow was measured T (°C)	Viscosity Correction factor F	Temperature in °C when permeate flow was measured T (°C)	Viscosity Correction factor F
4	1.567	25	0.89
5	1.519	26	0.871
6	1.472	27	0.851
7	1.428	28	0.833
8	1.386	29	0.815
9	1.346	30	0.798
10	1.307	31	0.781
11	1.271	32	0.765
12	1.235	33	0.749
13	1.202	34	0.734
14	1.169	35	0.719
15	1.139	36	0.705
16	1.109	37	0.692
17	1.081	38	0.678
18	1.053	39	0.665
19	1.027	40	0.653
20	1	41	0.641
21	0.978	42	0.629
22	0.955	43	0.618
23	0.933	44	0.607
24	0.911	45	0.596

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
## 4. Integrity Test

Ensure that the system is thoroughly cleaned and the membrane is thoroughly wetted out. Drain the system as thoroughly as possible for integrity test. This document provides two test methods, one of which may be applied according to actual operating conditions.

### 1. Diffusion Flow Method

1. Close the permeate valve and retentate valve, and immerse both lines under water.
2. Pressurize the retentate end with compressed air to increase the inlet pressure of the system to 15psi, then close the vent valve to maintain a positive pressure in the hollow fiber membrane.
3. Monitor the change of pressure with time. If the pressure drop is greater than 1psi/min, the system pipeline may have leakage. Check and tighten the pipeline, re-pressurize to 15psi, and monitor the pressure drop again.
4. After confirming that the integrity of the system is good, open the permeate valve, discharge the excess water, re-pressurize to 15psi, and measure the gas flow rate at the permeate end with a graduated cylinder full of water.
5. If the diffusion flow rate is less than 20ml/min/m<sup>2</sup>, it is considered that the integrity of the module is good. Release the pressure.
6. If the diffusion flow rate is greater than 20ml/min/m<sup>2</sup>, release the pressure, thoroughly wet the module again, and repeat step 1-5.


7. If the diffusion flow rate is still greater than 20ml/min/m<sup>2</sup>, it is considered that the membrane may have integrity defects. Do not use it.

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

## 2. Pressure Retaining Method

1. Close the permeate valve and retentate valve.
2. Slowly pump air by a peristaltic pump at the inlet end to increase the inlet pressure of the system to 15psi, then close the inlet valve to maintain the positive pressure in the hollow fiber membrane.
3. Monitor the change of pressure with time. If the pressure drop is greater than 1psi / min, the system pipeline may have leakage. Check and tighten the pipeline, re-pressurize to 15psi, and monitor the pressure drop again.
4. After confirming that the integrity of the system is good, open the retentate valve and permeate valve to drain the air and water in the system.
5. Close the retentate valve and keep the permeate valve open. Slowly pump air by the peristaltic pump at the inlet end to increase the inlet pressure of the system to 15psi, then close the liquid inlet valve and monitor the inlet pressure drop rate.
6. If the pressure drop is less than 2psi / min, it is considered that the integrity of the module is good for use. Release the pressure.
7. If the pressure drop is less than 2psi / min, release the pressure, fully wet out the components again, and repeat steps 1-5.
8. If the pressure drop is less than 2psi / min, it is considered that the membrane may have integrity defects. Do not use it.

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
## 5. Preprocessing Buffer Equilibrium (Optional)

The purpose of flushing hollow fiber membrane modules with buffer is to ensure that the pH and salt ion conditions of the system are as similar as possible to the initial state of the sample to be treated, so as to reduce membrane pollution and product loss.


1. Open both the retentate valve and the permeate valve; place both the retentate line and permeate line back to the feed tank.
2. Add sufficient buffer to the feed tank and circulate for 5-10min. Check pH or conductivity of the filtrate to ensure a similar condition to the buffer.

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
## 6. Cleaning After Use

The modules must be cleaned before storage and reuse. They should be flushed with buffer and/or DI water, and then washed with appropriate sanitizing agent. Under normal conditions, 0.5N NaOH may be used to circulate for 30min at 20-40°C. 300ppm NaCl may also be added to 0.5N NaOH solution if the cleaning effect is poor.


**Note:** Single use of Cobetter's hollow fiber cartridges is highly recommended. The cleaning method does not guarantee a full recovery of membrane performance. There might be compatibility issues if a specific solution/reagent is applied, or a special cleaning agent is needed for cleaning. Please consult Cobetter technical engineers for information before use.

1. Add sufficient buffer to the feed tank.
2. Close the permeate valve, open the retentate valve, and place the retentate and permeate line into the waste tank.
3. Start the pump slowly until the washing liquid fills the module, increase the pump speed and adjust the flow rate to 50% of the process flow rate.
4. A flushing volume of greater than 20L/m<sup>2</sup> is recommended.
5. Open the permeate valve, move the retentate line into the feed tank, adjust the retentate valve to let liquid flow out of the permeate end. The recommended permeate volume is greater than 20 L/m<sup>2</sup>. Drain the liquid in the system.
6. Add sufficient DI water to the cleaned feed tank and repeat step 2-5.
7. Add 20L/m<sup>2</sup> sanitizing agent into the cleaned feed tank.
8. Open the permeate valve and the retentate valve, and place the retentate and the permeate line into the waste tank.
9. Slowly start the pump until the sanitizing agent fills the system, increase the pump speed and adjust the flow rate to 100-150% of the process flow, so that 10% volume of sanitizing agent passes through the retentate and the permeate line, respectively.

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10. Move the retentate and permeate line into the feed tank, adjust the retentate valve so that there is liquid flowing out of the permeate end. Clean for 30min.
11. Empty the system and stop the pump.
12. Add sufficient DI water to the cleaned feed tank, rinse the retentate and permeate line until the pH of the effluent is neutral or the conductivity is low enough.

## 7. Hollow Fiber Storage

The hollow fiber module may be stored in 0.1N NaOH solution. For long-term storage (1-12 months), it is recommended to be stored in a refrigerated state (4°C).

Note: the appropriate storage of the modules ensures a wet membrane state and a effective prevention of biological pollution. Change of the storage solution once a month is highly recommended.

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