

Lifecube SAP Single-Use PUPSIT Assembly

Lifecube SAP system products are designed to perform PUPSIT (Pre-Use and Post-Sterilization Integrity Test) for filters in sterile filtration process and can be applied for single or two filters/redundant sterile filtration.



Regulatory Compliance: EU GMP Annex 1

PUPSIT can verify the integrity of a sterilizing grade filter before use and after sterilization.

Full Flexibility of Design, Limited Product Loss

Full line of single-use products including sterilizing-grade filters, storage bags, tubing, aseptic connectors, fittings, pressure gauge, in-line pressure sensors and monitor, etc. Meet almost all need for sterile filtration in final filling.

Delivery Time

Pre and self-produced assemblies, save at least 50% waiting times.

Standardized and Customized Solutions

Configuration-to-Order (CTO): Short delivery time
Engineer-To-Order (ETO) : Flexible design
Meet special single-use assemblies needs

In addition

- The design fully considers the process conditions and facility layout, such as WFI washing, pressure, air/water source etc. We will recommend the bag volume, tube size or combination solutions per the above factors
- Perform PUPSIT for filter with high bubble point
- Flushable components are pre-rinsed, effectively reduce the level of particulate matter and extracts
- Double-layer packaging, pre-irradiated, special needs can be customized three layer packaging, to meet customer needs in different clean areas of unpacking
- Self-produced components to ensure good quality traceability
- Various supporting documents, including validation guide, extractable data, system design guide, user manual, application data, etc


Quality assurance

- 100% Integrity test
- ISO 9001:2015 Quality Management System
- ISO Class 4.8 clean area assemble (equivalent to GMP Class A cleanliness)
- USP<87> Biological Reactivity Test, In Vitro
- USP<88> Class VI plastics Biological Reactivity Test, In Vivo
- Bacterial endotoxin meets WFI requirements (< 0.25 EU/mL)
- Particles matter meet the requirements in CP and USP<788> for large-volume parenterals
- Gamma radiation dose validation according to ISO11137, sterile packaging


Ordering Information

Customized PUPSIT designs

Single sterile Filtration
PUPSIT designs



Accessory
Stainless steel Holder




Key performance validation

- Leak Test
- Pressure Resistance Test
- Connect strength test


Typical Application

- Terminal Sterilization Filtration (EU GMP requirements)

Redundant sterile Filtration
PUPSIT designs (A)



Redundant sterile Filtration (B)



cobetter®
— filtration —

Tel 400-070-4266
Fax 0571-87704256
zh.cobetterfiltration.com

Please contact us for more information
Hangzhou Cobetter Filtration Equipment Co.,Ltd.

Sales Add 13-18 F, Building 1, Zicheng International Innovation Center, No.39 Jincheng Road,
Xiaoshan District, Hangzhou 311215, China

Factory Cobetter Park, Heshang New Material Industrial Park,
Xiaoshan District, Hangzhou 311265, China

Specifications subject to change without notice. CBPLSUSPA245b E